

BL40XU: Microbeam small-angle X-ray diffraction of hair

Noboru Ohta

1. Learning about BL40XU

- A) Light source, slits and mirror systems
- B) Microbeam small-angle X-ray diffraction

2. Training on adjustment of components and timing


- A) adjust the positions of pinholes and a beam stop
- B) adjust the position of a microscope for sample positioning
- C) find proper timing of an x-ray shutter and detector
- D) measure a beam size

3. Microbeam small-angle X-ray diffraction measurement

- A) estimate the pixel size of a 2D detector
- B) data acquisition of calibration sample
- C) scanning measurements using samples such as hair

4. Practicing data analysis

- A) handling of scattering data in a 2D array
 - read and display
 - slicing
 - circular averaging (partially)
 - find the beam center
- B) calibrate the scattering vector
- C) sequential processing



Schematic diagram of the microbeam small-angle x-ray diffraction set-up at BL40XU.

1. synchrotron x-ray source, 2. two focusing mirrors, 3. x-ray shutter, 4. collimating pinhole, 5. guard pinhole, 6. sample, 7. optical mirror for sample positioning, 8. microscope system, 9. vacuum pipe, 10. beam stop, 11. x-ray detector